




**FASTEST/MOST EFFECTIVE WAY TO COMBAT SKIN INFECTIONS!**

## PROVIDE THE SAFEST POSSIBLE ENVIRONMENT FOR COMPETITION

THE STERILASER DELIVERS 57,000 MICROWATTS OF GERMICIDAL ULTRAVIOLET (UV-C) LIGHT KILLING 99.5% OF ALL HARMFUL BACTERIA, FUNGI, AND VIRUSES IN LESS THAN 0.25 SECONDS. ALL COMMON PATHOGENS FOUND ON MATS THAT CAUSE DISEASES SUCH AS RINGWORM, IMPETIGO, AND MRSA ARE KILLED IN ONE PASS OF THE STERILASER. WHEN PULLED AT A NORMAL WALKING PACE, IT COMPLETELY STERILIZES AN ENTIRE 40' X 40' WRESTLING MAT IN 4 MINUTES



### KEY FEATURES

4 POWERFUL GERMICIDAL LAMPS



48" WIDTH/OPTIONAL BATTERY PACK



AUTOMATIC SHUT OFF IF TILTED



SHATTERPROOF COATED



# Compare for Yourself



## STERILASER

4 Minutes

Can be used between matches or rounds - as frequently as needed or desired.

Can be used by anyone walking at a normal pace

Eliminates 99.7% of pathogens with one pass

Easily rolls into storage closet

**Time to Sterilize Entire Mat**

**Frequency of Use**

**Knowledge Level**

**Effectiveness**

**Storage**



## MOP & BUCKET

45 Minutes

Can only be applied once before a match or tournament.

Proper use includes knowledge of the solution, application and disposal. insufficient knowledge of the process makes the application useless.

Virtually ineffective if not applied correctly.

Dumping mat cleaner down the drain is in violation of Federal EPA law.